

1 Collect your sample

DUST SAMPLER



Wipe a visibly dusty surface with the absorbent pad of the CardSampler. Keep the surface dry - do not wet it while collecting.

2 Return the sampler for lab analysis

- Seal the sampler in a ziplock bag to protect it during transit.
- Fill in the submission form below, then cut it out and place it in the envelope with the sampler.
- Cut out the return shipping label and tape it to the envelope - cover the whole label with clear tape so it stays secure.

Once your sampler arrives, we will email you a quote for the analysis.

CUT OUT THE TWO BOXES BELOW

RETURN SHIPPING LABEL

CUSTOMS DECLARATION

Harmless material sample
with no commercial value.

RETURN TO

Technology Care GmbH
Birmensdorferstrasse 467
8055 Zurich
SWITZERLAND

SUBMISSION FORM

Your Email

Sampler ID

Your Name & Address

GOOD TO KNOW

- Handle the sampler by its edges only - avoid touching the collection surface.
- Keep the start date legible; it sets the exposure period for your results.
- Questions? Visit www.technologycare.com or contact us before returning your sampler.